

MAGIQ Sphere

MRI Phantoms

Features

- Spherical flood-fill phantom
- Choice of internal features; no features (hollow), test tube holder or features as per MAGIQ
- Sized to fit head coils; outer diameter 190mm
- Phantom support

Benefits

- Multiple plane testing for Uniformity and SNR without repositioning of the phantom
- Spherical shape is less susceptible to Gibbs artifact (a series of lines in the MR image parallel to abrupt and intense changes in the signal) than cylindrical phantoms
- Simple support prevents accidental rolling of phantom when rested on patient couch



MAGIQ Sphere (hollow, no features) positioned on phantom support