



VISIRAY

Radiopaque & Phosphorescent Rulers

Radiopaque and phosphorescent rulers for quality control of congruence between light, radiation and x-ray fields.

Features:

- Radiopaque markings (1mm)
 - Visible markings (1mm)
 - Phosphor coating (28-120 kV)
- Removable protective casing

Specification:

- Size: 45cm, 30cm, 20cm and 10cm
- Bespoke sizes on request





VISIRAY

Radiopaque & Phosphorescent Rulers




Protective cases can be removed by removing the magnetically attached cover and extracting the ruler.



VISIRAY

Radiopaque & Phosphorescent Rulers



Frame captured from 4K resolution camera
showing radiation field size




**4K Camera
Available**



VISIRAY

Radiopaque & Phosphorescent Rulers



Radiograph showing x-ray image field size